

Corning® Aspirator

A convenient, user-friendly device that facilitates the aspiration of liquids from lab containers

CORNING

The Corning Aspirator is an aspiration tool used for the removal of liquids from a variety of lab vessels. The device allows precise fingertip aspiration control, and is a sensible alternative to the use of glass pipets for supernatant disposal.

Designed for a comfortable fit in either hand, the aspirator enables easy handling of liquids from a variety of collection devices such as tubes, dishes, bottles, T-flasks, and multiwell plates.

- ▶ Vacuum regulation controlled by a touch-sensitive button allows for quick and comfortable aspiration of liquids
- ▶ Ergonomic design makes the device perfect for long work sessions or repetitive aspiration tasks
- ▶ The aspirator comes with a rubber adapter for accessory attachment, as well as a single-channel adapter for 1 to 200 µL disposable pipet tips
- ▶ An 8-channel adapter is available for working with 96-well microplates
- ▶ Anti-drip technology prevents liquids leaking from the pipet tips, reducing contamination risks and the likelihood of collecting erroneous results
- ▶ Entirely autoclavable for maximum safety when working with hazardous materials

Ordering Information

Products may not be available in all markets.

Cat. No.	Description	Qty/Cs
4930	Corning Aspirator	1
155850	Single-channel adapter for the Corning Aspirator	1
155505	Rubber adapter for glass pipets	1

Spare Parts and Accessories (sold separately)


Cat. No.	Description	Qty/Cs
4931	8-channel adapter for the Corning Aspirator	1
155605	Tappet grey with O-ring, valve plunger	1
151230	Plunger spring	1
155805	O-ring, 2.57 x 1.78 mm FKM 60 shore A for valve plunger	1
155800	O-ring, 9.25 x 1.78 mm FKM for base screw	1
155830	O-ring, 6.75 x 1.78 mm FKM for end piece	1
155855	Quick-locking screw, orange	1



8-channel adapter

Related Products

Products may not be available in all markets.



Cat. No.	Description	Qty/Pk	Qty/Cs
4860	Corning® 1 - 200 µL universal fit racked pipet tips, yellow, polypropylene, sterile	96	960
4121	Corning DeckWorks 1 - 200 µL pipet tips, graduated, hinged racks, natural, polypropylene, sterile	960	3,840
4144	Corning DeckWorks 1 - 200 µL pipet tip station, graduated, natural, polypropylene, sterile	480	4,800
4020	Stripette® 2 mL serological pipets, polystyrene, bulk packed, sterile	50	1,000
4486	Stripette 2 mL serological pipets, polystyrene, individually paper/plastic wrapped, sterile	50	1,000
4021	Stripette 2 mL serological pipets, polystyrene, individually plastic wrapped, sterile	100	1,000
9186	Costar® 2 mL aspirating pipets, polystyrene, without graduations, individually wrapped, sterile	1	1,000

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated
Life Sciences

www.corning.com/lifesciences

NORTH AMERICA

t 800.492.1110
t 978.442.2200

ASIA/PACIFIC

Australia/New Zealand

t 61 427286832

Chinese Mainland

t 86 21 3338 4338

India

t 91 124 4604000

Japan

t 81 3-3586 1996

Korea

t 82 2-796-9500

Singapore

t 65 6572-9740

Taiwan

t 886 2-2716-0338

EUROPE

CSEurope@corning.com

France

t 0800 916 882

Germany

t 0800 101 1153

The Netherlands

t 020 655 79 28

United Kingdom

t 0800 376 8660

All Other European Countries

t +31 (0) 206 59 60 51

LATIN AMERICA

grupoLA@corning.com

Brazil

t 55 (11) 3089-7400

Mexico

t (52-81) 8158-8400