




# Automation Tips

Overview – Automation Tips . . . . .	3
Aurora Biomed . . . . .	5
Agilent . . . . .	6
Agilent BioTek . . . . .	7
Beckman Coulter . . . . .	8
Corbett . . . . .	10
Dynamic Devices . . . . .	11
Hamilton . . . . .	13
Molecular Devices . . . . .	16
PerkinElmer . . . . .	17
Tecan . . . . .	19
Axygen® Automation Reservoirs. . . . .	23

## Product Ordering Information

Products may not be available in all markets. For information on Purchasing Options, Terms and Conditions of Sale, Return and Repair Policies, and Warranty/Guarantee Registration, visit [www.corning.com/how-to-buy](http://www.corning.com/how-to-buy).

# Overview – Automation Tips



Axygen® automation tips are manufactured to meet precise specifications required for accurate sample delivery in automated applications. Our tips are manufactured using high quality molds coupled with Corning process controls, ensuring complete compatibility and consistent fit. These tips are manufactured using automation to ensure tips are free from contamination and defects that may be caused by hand racking.

- ▶ Wide product offering combined with in-house engineering expertise guarantees optimum workstation performance
- ▶ Complete control of manufacturing process with comprehensive QA and QC testing processes ensures defects are eliminated
- ▶ Optimized fine point for accurate dispensing
- ▶ Filter barrier prevents aerosol and liquid contamination
- ▶ “Plug and Play” – no software or hardware changes needed
- ▶ Color-coded racks for easy identification of tip size
- ▶ ANSI/SLAS footprint rack dimensions for precise fit on the deck (except where noted)
- ▶ Nonsterile, sterile, filter/sterile versions are available
- ▶ RNase-/DNase-free, human gDNA-free
- ▶ Nonpyrogenic

## Axygen Automation Tips are Compatible with:

### Aurora Biomed

- ▶ VERSA™

### Agilent

- ▶ VPrep® and Bravo

### Agilent BioTek

- ▶ Precision™ and Precision XS

### Beckman Coulter

- ▶ Biomek® FX, NX, i5, i7

### Corbett

- ▶ CAS
- ▶ X-tractor Gene, QiaXtractor

### Dynamic Devices

- ▶ Oasis LM

### Hamamatsu

- ▶ 6000 and 7000 Series

### Hamilton

- ▶ CO-RE™ Microlab® Prep™, STAR™ Line, Nimbus®, and VANTAGE®

### Molecular Devices

- ▶ PatchXpress, FLIPR® Tetra, and FlexStation®

### PerkinElmer

- ▶ Janus® and Evolution P3
- ▶ VariSpan

### Qiagen


- ▶ BioRobot®

### Tecan

- ▶ MCA96 Head
- ▶ LiHa Head

### 5 Reasons to Choose Axygen® Automation Tips


- 1 Automated Assembly**  
Ensures product is free from contamination.
- 2 Electronic Monitoring of Manufacturing**  
Advanced technology supported by experienced personnel results in high quality products.
- 3 In-house Testing**  
In-house automation lab allows product to be validated prior to release.
- 4 Customer Site Testing of New Products**  
All new products are validated at customer locations to ensure proper function.
- 5 Broad Product Offering**  
Corning routinely offers automation tips not offered by the robotic workstation manufacturer.



## Axygen® 96-well Format Tips

Compatible with Aurora Biomed VERSA™ Workstations

- ▶ Compatible head: 8-channel
- ▶ Non-ANSI/SLAS hanging tip, polycarbonate rack



### Clear

Cat. No.	Description	Packaging
TTF-20-C-HTR-S	10 µL clear filtered tips, sterile	96 tips/rack, 24 racks/case
TTF-50-C-HTR-S	25 µL clear filtered tips, sterile	96 tips/rack, 24 racks/case
TT-50-C-HTR	50 µL clear non-filtered tips	96 tips/rack, 24 racks/case
TT-50-C-HTR-36	50 µL clear non-filtered tips, bulk pack	96 tips/rack, 2 racks/pk, 36 pk/case
TTF-200-C-HTR-S	175 µL clear filtered tips, sterile	96 tips/rack, 24 racks/case
TT-200-C-HTR	200 µL clear non-filtered tips	96 tips/rack, 24 racks/case
TT-200-C-HTR-S	200 µL clear non-filtered tips, sterile	96 tips/rack, 24 racks/case
TT-1000-C-HTR	1000 µL clear non-filtered tips	96 tips/rack, 16 racks/case
TT-1000-C-HTR-S	1000 µL clear non-filtered tips, sterile	96 tips/rack, 16 racks/case
TTF-1000-C-HTR-S	1000 µL clear filtered tips, sterile	96 tips/rack, 16 racks/case



### Conductive Black

Cat. No.	Description	Packaging
TTF-20-CBK-HTR-S	10 µL black filtered tips, sterile	96 tips/rack, 24 racks/case
TT-20-CBK-HTR	20 µL black non-filtered tips	96 tips/rack, 24 racks/case
TTF-50-CBK-HTR-S	25 µL black filtered tips, sterile	96 tips/rack, 24 racks/case
TT-50-CBK-HTR	50 µL black non-filtered tips	96 tips/rack, 24 racks/case
TT-50-CBK-HTR-36	50 µL black non-filtered tips, bulk pack	96 tips/rack, 2 racks/pk, 36 pk/case
TT-50-CBK-HTR-S	50 µL black non-filtered tips, sterile	96 tips/rack, 24 racks/case
TTF-200-CBK-HTR-S	175 µL black filtered tips, sterile	96 tips/rack, 24 racks/case
TT-200-CBK-HTR	200 µL black non-filtered tips	96 tips/rack, 24 racks/case
TT-200-CBK-HTR-S	200 µL black non-filtered tips, sterile	96 tips/rack, 24 racks/case
TT-1000-CBK-HTR	1000 µL black non-filtered tips	96 tips/rack, 16 racks/case
TT-1000-CBK-HTR-24	1000 µL black non-filtered tips, bulk pack	96 tips/rack, 2 racks/pk, 24 pk/case
TT-1000-CBK-HTR-S	1000 µL black non-filtered tips, sterile	96 tips/rack, 16 racks/case
TTF-1000-CBK-HTR-S	1000 µL black filtered tips, sterile	96 tips/rack, 16 racks/case

### Axygen® 96-well Format Tips, Clear

Compatible with Agilent VPrep® and Bravo Workstations

- ▶ Compatible heads: 96LT, 8LT
- ▶ ANSI/SLAS, polypropylene rack




Cat. No.	Description	Packaging
VT-20-R	20 µL clear non-filtered tips	96 tips/rack, 50 racks/case
VTF-20-R-S	20 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
VT-50-R	50 µL clear non-filtered tips	96 tips/rack, 50 racks/case
VTF-165-R-S	165 µL clear filtered tips, sterile (Old V-11 Style)	96 tips/rack, 50 racks/case
VTF-205-R-S	200 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
VT-250-R	250 µL clear non-filtered tips (Old V-11 Style)	96 tips/rack, 50 racks/case
VT-250-R-S	250 µL clear non-filtered tips, sterile (Old V-11 Style)	96 tips/rack, 50 racks/case
VT-255-R	250 µL clear non-filtered tips	96 tips/rack, 50 racks/case
VT-255-R-S	250 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case

### Axygen 384-well Format Tips, Clear

Compatible with Agilent VPrep and Bravo Workstations

- ▶ Compatible heads: 384ST, 96ST, 16ST
- ▶ ANSI/SLAS, conductive black polypropylene rack






Cat. No.	Description	Packaging
VT-384-5UL-R	5 µL clear non-filtered tips	384 tips/rack, 50 racks/case
VT-384-10UL-R	10 µL clear non-filtered tips	384 tips/rack, 50 racks/case
VT-384-10UL-R-S	10 µL clear non-filtered tips, sterile	384 tips/rack, 50 racks/case
VTF-384-15UL-R-S	15 µL clear filtered tips, sterile	384 tips/rack, 50 racks/case
VT-384-31UL-R	30 µL clear non-filtered tips	384 tips/rack, 50 racks/case
VT-384-31UL-R-S	30 µL clear non-filtered tips, sterile	384 tips/rack, 50 racks/case
VTF-384-50UL-R-S	50 µL clear filtered tips, sterile	384 tips/rack, 50 racks/case
VT-384-70UL-R	70 µL clear non-filtered tips	384 tips/rack, 50 racks/case
VT-384-70UL-R-S	70 µL clear non-filtered tips, sterile	384 tips/rack, 50 racks/case

## Axygen® 96-well Format Tips, Clear

Compatible with BioTek Precision™ and Precision XS Workstations

- ▶ Compatible heads: 8- or 12-channel
- ▶ ANSI/SLAS, polypropylene rack




Cat. No.	Description	Packaging
EV-50-R	50 µL clear non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
EV-50-R-S	50 µL clear non-filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
EVF-50-R-S	50 µL clear filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
EVF-100-R-S	80 µL clear filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
EV-100-R	100 µL clear non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
EVF-180-R-S	180 µL clear filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
EVF-205-WB-R-PB	180 µL clear wide bore filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
EVF-205-WB-R-S	180 µL clear wide bore filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
EV-200-R	200 µL clear non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
EV-200-R-S	200 µL clear non-filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case

### Axygen® 96-well Format Tips, Clear

Compatible with Beckman Biomek® FX, NX, i5 and i7 Workstations

- ▶ Compatible heads: AP-96, MP20, MP200
- ▶ ANSI/SLAS, polypropylene rack




Cat. No.	Description	Packaging
FX-10-R	10 µL clear non-filtered tips	96 tips/rack, 50 racks/case
FX-20-R	20 µL clear non-filtered tips	96 tips/rack, 50 racks/case
FX-20-R-S	20 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case
FXF-20-R-S	20 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
FXF-50-R-S	30 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
FX-50-R	50 µL clear non-filtered tips	96 tips/rack, 50 racks/case
FX-50-R-S	50 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case
FXF-180-R-S	165 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
FXF-165-WB-R-S*	165 µL clear wide bore filtered tips, sterile	96 tips/rack, 50 racks/case
FXF-200-R-S	200 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
FXF-205-WB-R-S	200 µL clear wide bore filtered tips, sterile	96 tips/rack, 50 racks/case
FX-250-R	250 µL clear non-filtered tips	96 tips/rack, 50 racks/case
FX-250-R-S	250 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case
FX-255-WB-R	250 µL clear wide bore non-filtered tips	96 tips/rack, 50 racks/case

\*Made to order. Please contact your local Corning Account Manager for more details.

### Axygen® 384-well Format Tips, Clear

Compatible with Beckman Biomek® FX, NX, i5 and i7 Workstations

- ▶ Compatible heads: AP-384, 384-channel low volume (FXF-384-XL-R-S only)
- ▶ ANSI/SLAS, conductive black polypropylene rack




Cat. No.	Description	Packaging
FX-384-R	30 µL clear non-filtered tips	384 tips/rack, 50 racks/case
FX-384-R-S	30 µL clear non-filtered tips, sterile	384 tips/rack, 50 racks/case
FX-384-30FP-R	30 µL clear fine point non-filtered tips	384 tips/rack, 50 racks/case
FX-384-30FP-R-S	30 µL clear fine point non-filtered tips, sterile	384 tips/rack, 50 racks/case
FXF-384-30FP-R-S	30 µL clear fine point filtered tips, sterile	384 tips/rack, 50 racks/case
FXF-384-XL-R-S	30 µL clear extended length filtered tips, sterile	384 tips/rack, 50 racks/case
FX-384-45FP-R	45 µL clear fine point non-filtered tips	384 tips/rack, 50 racks/case
FXF-384-45FP-R-S	45 µL clear fine point filtered tips, sterile	384 tips/rack, 50 racks/case
FX-384-XL-R	50 µL clear extended length non-filtered tips	384 tips/rack, 50 racks/case
FX-384-XL-R-S	50 µL clear extended length non-filtered tips, sterile	384 tips/rack, 50 racks/case




### Axygen® 96-well Format Tips

Compatible with Corbett CAS Workstations

- ▶ Compatible heads: CAS-1200 single-channel (20 - 150 µL/sec.)
- ▶ ANSI/SLAS hanging tip, polycarbonate rack



#### Clear



Cat. No.	Description	Packaging
TT-50-C-HTR	50 µL clear non-filtered tips	96 tips/rack, 24 racks/case
TT-50-C-HTR-36	50 µL clear non-filtered tips, bulk pack	96 tips/rack, 2 racks/pk, 36 pk/case
TTF-50-C-HTR-S	25 µL clear filtered tips, sterile	96 tips/rack, 24 racks/case
TTF-200-C-HTR-S	175 µL clear filtered tips, sterile	96 tips/rack, 24 racks/case
TT-200-C-HTR	200 µL clear non-filtered tips	96 tips/rack, 24 racks/case
TT-200-C-HTR-S	200 µL clear non-filtered tips, sterile	96 tips/rack, 24 racks/case

#### Conductive Black



Cat. No.	Description	Packaging
TTF-50-CBK-HTR-S	25 µL black filtered tips, sterile	96 tips/rack, 24 racks/case
TT-50-CBK-HTR	50 µL black non-filtered tips	96 tips/rack, 24 racks/case
TT-50-CBK-HTR-36	50 µL black non-filtered tips, bulk pack	96 tips/rack, 2 racks/pk, 36 pk/case
TT-50-CBK-HTR-S	50 µL black non-filtered tips, sterile	96 tips/rack, 24 racks/case
TT-200-CBK-HTR	200 µL black non-filtered tips	96 tips/rack, 24 racks/case
TT-200-CBK-HTR-S	200 µL black non-filtered tips, sterile	96 tips/rack, 24 racks/case
TTF-200-CBK-HTR-S	200 µL black filtered tips, sterile	96 tips/rack, 24 racks/case

### Axygen 96-well Format Tips, Clear

Compatible with Corbett X-tractor Gene and QiaXtractor Workstations

- ▶ Compatible heads: X-tractor Gene 8-channel, CAS-4200, CAS-4800, QiaXtractor
- ▶ Standard ANSI/SLAS, polypropylene (CTF) rack




Cat. No.	Description	Packaging
CTF-199-R-S	200 µL clear filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case

### Axygen® 96-well Format Tips, Clear

Compatible with Dynamic Devices Oasis LM Workstations

- ▶ Compatible heads: 8- or 96-channel standard volume
- ▶ ANSI/SLAS, polypropylene rack




Cat. No.	Description	Packaging
EV-50-R	50 µL clear non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
EV-50-R-S	50 µL clear non-filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
EVF-50-R-S	50 µL clear filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
EVF-100-R-S	80 µL clear filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
EV-100-R	100 µL clear non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
EVF-180-R-S	180 µL clear filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
EVF-205-WB-R-PB	180 µL clear wide bore filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
EVF-205-WB-R-S	180 µL clear wide bore filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
EV-200-R	200 µL clear non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
EV-200-R-S	200 µL clear non-filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case

### Axygen 96-well Format Tips, Clear, Paper Bottom

Compatible with Dynamic Devices Oasis LM Workstations

- ▶ Compatible heads: 8- or 96-channel standard volume
- ▶ ANSI/SLAS, polypropylene rack




Cat. No.	Description	Packaging
OTF-50-R-PB	50 µL clear filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
OTF-100-R	100 µL clear non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
OTF-100-R-PB	100 µL clear filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
OTF-180-R-PB	180 µL clear filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
OTF-200-R	200 µL clear non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case

## Dynamic Devices (continued)

### Axygen® 384-well Format Tips, Clear

Compatible with Dynamic Devices Oasis LM Workstations

- Compatible heads: 384-channel low volume
- ANSI/SLAS, polypropylene rack




Cat. No.	Description	Packaging
FX-384-R	30 µL clear non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
FX-384-R-S	30 µL clear non-filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
FX-384-30FP-R	30 µL clear fine point non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
FX-384-30FP-R-S	30 µL clear fine point non-filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
FXF-384-30FP-R-S	30 µL clear fine point filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
FXF-384-XL-R-S	30 µL clear filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
FX-384-45FP-R	45 µL clear fine point non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
FXF-384-45FP-R-S	45 µL clear fine point filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case
FX-384-XL-R	50 µL clear non-filtered tips	96 tips/rack, 10 racks/pk, 5 pk/case
FX-384-XL-R-S	50 µL clear non-filtered tips, sterile	96 tips/rack, 10 racks/pk, 5 pk/case

## Hamamatsu

### Axygen 384-well Format Tips, Black Conductive

Compatible with Hamamatsu 6000 and 7000 Series Workstations

- Compatible heads: FDSS 6000, FDSS 7000
- Non-ANSI/SLAS hanging tip, polypropylene rack





Cat. No.	Description	Head Compatibility	Packaging
FDSS-384-BK-R*	50 µL black non-filtered tips	FDSS 6000	96 tips/rack, 10 racks/pk, 5 pk/case
FDSS-1536-BK-R*	50 µL black non-filtered tips	FDSS 7000	96 tips/rack, 10 racks/pk, 5 pk/case

\*Made to order. Please contact your local Corning Account Manager for more details.

## Axygen® 96-well Format Tips, Hanging Tip Rack

Compatible with Hamilton CO-RE™ Microlab® Prep™, STAR™ Line, Nimbus®, and VANTAGE® Workstations

- ▶ Compatible heads: 8- and 96-channel
- ▶ ANSI/SLAS, hanging tip, polypropylene rack



### Conductive Black

Cat. No.	Description	Qty/Pk	Qty/Cs
HT-10-CBK-HTR	10 µL liquid level sensing automation tips	96	2,304
HT-10-CBK-HTR-5	10 µL liquid level sensing automation tips	480	5,760
HT-10-CBK-HTR-S-5	10 µL liquid level sensing automation tips, sterile	480	5,760
HTF-10-CBK-HTR-S	10 µL liquid level sensing filtered automation tips, sterile	96	2,304
HTF-10-CBK-HTR-S-5	10 µL liquid level sensing filtered automation tips, sterile	480	5,760
HT-50-CBK-HTR	50 µL liquid level sensing automation tips	96	2,304
HT-50-CBK-HTR-5	50 µL liquid level sensing automation tips	480	5,760
HT-50-CBK-HTR-S-5	50 µL liquid level sensing automation tips, sterile	480	5,760
HTF-50-CBK-HTR-S	50 µL liquid level sensing filtered automation tips, sterile	96	2,304
HTF-50-CBK-HTR-S-5	50 µL liquid level sensing filtered automation tips, sterile	480	5,760
HT-300-CBK-HTR	300 µL liquid level sensing automation tips	96	2,304
HT-300-CBK-HTR-5	300 µL liquid level sensing automation tips	480	5,760
HT-300-CBK-HTR-S	300 µL liquid level sensing automation tips, sterile	96	2,304
HT-300-CBK-HTR-S-5	300 µL liquid level sensing automation tips, sterile	480	5,760
HTF-300-CBK-HTR-S	300 µL liquid level sensing filtered automation tips, sterile	96	2,304
HTF-300-CBK-HTR-S5	300 µL liquid level sensing filtered automation tips, sterile	480	5,760
HT-1000-CBK-HTR-5	1000 µL liquid level sensing automation tips	480	3,840
HT-1000-CBK-HTR-S	1000 µL liquid level sensing automation tips, sterile	96	1,536
HT-1000-CBK-HTR-S5	1000 µL liquid level sensing automation tips, sterile	480	3,840
HTF-1000-CBK-HTR-5	1000 µL liquid level sensing filtered automation tips	480	3,840
HTF-1000-CBK-HTR-S	1000 µL liquid level sensing filtered automation tips, sterile	96	1,536
HTF-1000-CBK-HTR-S5	1000 µL liquid level sensing filtered automation tips, sterile	480	3,840



**Axygen® 96-well Format Tips, Hanging Tip Rack (continued)**

Compatible with Hamilton CO-RE™ Microlab® Prep™, STAR™ Line, Nimbus®, and VANTAGE® Workstations (continued)


**Clear**

Cat. No.	Description	Qty/Pk	Qty/Cs
HT-10-C-HTR-5	10 µL automation tips	480	5,760
HT-10-C-HTR-S-5	10 µL automation tips, sterile	480	5,760
HTF-10-C-HTR-S-5	10 µL filtered automation tips, sterile	480	5,760
HT-50-C-HTR-5	50 µL automation tips	480	5,760
HT-50-C-HTR-S-5	50 µL automation tips, sterile	480	5,760
HTF-50-C-HTR-S-5	50 µL filtered automation tips, sterile	480	5,760
HT-300-C-HTR-5	300 µL automation tips	480	5,760
HT-300-C-HTR-S-5	300 µL automation tips, sterile	480	5,760
HTF-300-C-HTR-S	300 µL filtered automation tips, sterile	96	2,304
HTF-300-C-HTR-S-5	300 µL filtered automation tips, sterile	480	5,760
HT-1000-C-HTR-5	1000 µL automation tips	480	3,840
HT-1000-C-HTR-S-5	1000 µL automation tips, sterile	480	3,840
HTF-1000-C-HTR-S	1000 µL filtered automation tips, sterile	96	1,536
HTF-1000-C-HTR-S-5	1000 µL filtered automation tips, sterile	480	3,840

**Axygen 96-well Format Tips, Nested Tip Rack**

Compatible with Hamilton CO-RE Microlab Prep, STAR Line, Nimbus, and VANTAGE Workstations

- ▶ Compatible heads: 8- and 96-channel
- ▶ ANSI/SLAS nested, polypropylene rack

**Conductive Black**

Cat. No	Description	Qty/Pk	Qty/Cs
HT-10-CBK-STK	10 µL liquid level sensing automation tips	384	11,520
HT-10-CBK-STK-S	10 µL liquid level sensing automation tips, sterile	384	11,520
HT-50-CBK-STK	50 µL liquid level sensing automation tips	384	11,520
HT-50-CBK-STK-S	50 µL liquid level sensing automation tips, sterile	384	11,520
HT-300-CBK-STK	300 µL liquid level sensing automation tips	384	11,520
HT-300-CBK-STK-S	300 µL liquid level sensing automation tips, sterile	384	11,520


**Clear**

Cat. No	Description	Qty/Pk	Qty/Cs
HT-10-C-STK	10 µL automation tips	384	11,520
HT-10-C-STK-S	10 µL automation tips, sterile	384	11,520
HT-50-C-STK	50 µL automation tips	384	11,520
HT-50-C-STK-S	50 µL automation tips, sterile	384	11,520
HT-300-C-STK	300 µL automation tips	384	11,520
HT-300-C-STK-S	300 µL automation tips, sterile	384	11,520

### Axygen® 384-well Format Clear Polypropylene Tips

Compatible with Hamilton CO-RE™ Microlab® Prep™, STAR™ Line, Nimbus®, and VANTAGE® Workstations

- ▶ Compatible head: 384-channel
- ▶ ANSI/SLAS, anti-static hanging tip, black conductive polypropylene rack



Cat. No	Description	Qty/Pk	Qty/Cs
FX-384-R	30 µL clear tips	3,840	19,200
FX-384-R-S	30 µL clear tips, sterile	3,840	19,200
FX-384-30FP-R	30 µL clear fine point tips	3,840	19,200
FX-384-30FP-R-S	30 µL clear fine point tips, sterile	3,840	19,200
FXF-384-30FP-R-S	30 µL clear fine point tips, filtered, sterile	3,840	19,200
FX-384-45FP-R	45 µL clear fine point tips	3,840	19,200
FX-384-45FP-R-S	45 µL clear fine point tips, sterile	3,840	19,200
FXF-384-45FP-R-S	45 µL clear fine point tips, filtered, sterile	3,840	19,200
FX-384-XL-R	50 µL clear extra long tips	3,840	19,200
FX-384-XL-R-S	50 µL clear extra long tips, sterile	3,840	19,200
FXF-384-XL-R-S	50 µL clear extra long filter tips, sterile	3,840	19,200

### Axygen® 96-well Format Tips, Clear and Black Conductive

Compatible with Molecular Devices PatchXpress, FLIPR® Tetra, and FlexStation® Workstations

- ▶ ANSI/SLAS (TT-50-C-HTR ANSI/SLAS hanging tips), polypropylene rack




Cat. No.	Description	Workstation Compatibility	Packaging
TT-50-C-HTR	50 µL black non-filtered tips	PatchXpress 7000A	96 tips/rack, 50 racks/case
PK-200-R	200 µL clear non-filtered tips	FLIPR Tetra	96 tips/rack, 50 racks/case
RMT-200-C-R	200 µL clear non-filtered tips	FlexStation	96 tips/rack, 50 racks/case
PK-200-BK-R	200 µL black non-filtered tips	FLIPR Tetra	96 tips/rack, 50 racks/case
RMT-200-BK-R	200 µL black non-filtered tips	FlexStation	96 tips/rack, 50 racks/case

### Axygen 384-well Format Tips, Clear and Black Conductive

Compatible with Molecular Devices FLIPR Tetra and FLIPR Liberty 384 Workstations

- ▶ ANSI/SLAS, polypropylene rack





Cat. No.	Description	Workstation Compatibility	Packaging
FT-384-C-R	30 µL clear non-filtered tips	FLIPR Tetra	384 tips/rack, 50 racks/case
FT-384-BK-R	30 µL black non-filtered tips	FLIPR Tetra	384 tips/rack, 50 racks/case
PK-384-R	30 µL clear non-filtered tips	FLIPR Tetra Type A	384 tips/rack, 50 racks/case
LT-384-R	30 µL clear non-filtered tips	FLIPR Liberty B (Liberty)	384 tips/rack, 50 racks/case

## Axygen® 96-well Format Clear MDT Tips

Compatible with PerkinElmer Janus® and Evolution P3 Workstations

- ▶ Compatible heads: MDT-96
- ▶ ANSI/SLAS, polypropylene rack



Cat. No.	Description	Packaging
PK-20-R	20 µL clear non-filtered tips	96 tips/rack, 50 racks/case
PK-50-R	50 µL clear non-filtered tips	96 tips/rack, 50 racks/case
PK-50-R-S	50 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case
PK-200-R	200 µL clear non-filtered tips	96 tips/rack, 50 racks/case

## Axygen 384-well Format Tips, Clear MDT

Compatible with PerkinElmer Janus and Evolution P3 Workstations

- ▶ Compatible heads: MDT-384-30 µL
- ▶ ANSI/SLAS, polypropylene rack




Cat. No.	Description	Packaging
PK-384-R	30 µL clear non-filtered tips	384 tips/rack, 50 racks/case
PK-384-R-S	30 µL clear non-filtered tips, sterile	384 tips/rack, 50 racks/case
PKF-384-R-S	30 µL clear filtered tips, sterile	384 tips/rack, 50 racks/case



### Axygen® 96-well Format Tips, Clear and Black Conductive

Compatible with PerkinElmer VariSpan Head

- ▶ ANSI/SLAS, polypropylene rack




Cat. No.	Description	Packaging
RRF-25-C-R-S	25 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
RRF-25-CBK-R-S	25 µL black liquid level sensing filtered tips, sterile	96 tips/rack, 50 racks/case
RRP-50-C-R	50 µL clear non-filtered tips	96 tips/rack, 50 racks/case
RRP-50-CBK-R	50 µL black liquid level sensing non-filtered tips	96 tips/rack, 50 racks/case
RRF-175-C-R-S	175 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
RRF-175-CBK-R-S	175 µL black liquid level sensing filtered tips, sterile	96 tips/rack, 50 racks/case
RRP-200-C-R	200 µL clear non-filtered tips	96 tips/rack, 50 racks/case
RRP-200-C-R-S	200 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case
RRP-200-CBK-R	200 µL black liquid level sensing non-filtered tips	96 tips/rack, 50 racks/case

### Qiagen

### Axygen 96-well Format Tips, Black Conductive

Compatible with Qiagen BioRobot® Series Workstations

- ▶ ANSI/SLAS hanging tip, polycarbonate rack




Cat. No.	Description	Packaging
Q-OCF-200-R-S	200 µL black filtered tips, sterile	96 tips/rack, 5 racks/case
QT-300-CBK-R	300 µL black non-filtered tips	96 tips/rack, 24 racks/case
QTF-1100-CBK-R-S	1000 µL black filtered tips, sterile	96 tips/rack, 16 racks/case

### Axygen® 96-well Format Tips, Clear

Compatible with Tecan MCA96 Head

▶ ANSI/SLAS, polypropylene rack




Cat. No.	Description	Packaging
EV-50-R	50 µL clear non-filtered tips	96 tips/rack, 50 racks/case
EV-50-R-S*	50 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case
EVF-50-R-S	50 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
EV-100-R	100 µL clear non-filtered tips	96 tips/rack, 50 racks/case
EV-100-R-S	100 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case
EVF-100-R-S*	100 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
EVF-180-R-S	180 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
EVF-205-WB-R-S	180 µL clear filtered tips, wide bore, sterile	96 tips/rack, 50 racks/case
EV-200-R	200 µL clear non-filtered tips	96 tips/rack, 50 racks/case
EV-200-R-S	200 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case

\*Made to order. Please contact your local Corning Account Manager for more details.

### Axygen 96-well Format Tips, Clear

Compatible with Tecan MCA96 Head (ANSI/SLAS Footprint Nested Rack)

▶ ANSI/SLAS nested, polypropylene rack




Cat. No.	Description	Packaging
EV-50-NTR	50 µL clear non-filtered tips	96 tips/rack, 8 racks/stack, 5 stacks/case
EV-50-NTR-S	50 µL clear non-filtered tips, sterile	96 tips/rack, 8 racks/stack, 5 stacks/case
EV-200-NTR	200 µL clear non-filtered tips	96 tips/rack, 8 racks/stack, 5 stacks/case
EV-200-NTR-S	200 µL clear non-filtered tips, sterile	96 tips/rack, 8 racks/stack, 5 stacks/case

## Axygen® 96-well Format Tips, Clear

Compatible with Tecan MCA96 Head

- ▶ Non-ANSI/SLAS SciClone® style, polypropylene rack




Cat. No.	Description	Packaging
ZT-50-R*	50 µL clear non-filtered tips	96 tips/rack, 50 racks/case
ZT-50-R-S*	50 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case
ZTF-50-R-S	50 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
ZT-100-R	100 µL clear non-filtered tips	96 tips/rack, 50 racks/case
ZT-100-R-S*	100 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case
ZTF-100-R-S	100 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
ZTF-200-R-S	180 µL clear filtered tips, sterile	96 tips/rack, 50 racks/case
ZT-200-R	200 µL clear non-filtered tips	96 tips/rack, 50 racks/case
ZT-200-R-S	200 µL clear non-filtered tips, sterile	96 tips/rack, 50 racks/case











\*Made to order. Please contact your local Corning Account Manager for more details.

### Axygen® 96-well Format Tips, Clear

Compatible with Tecan LiHa Head

- ▶ ANSI/SLAS hanging tip, polypropylene rack




	Cat. No.	Description	Packaging
	TTF-20-C-HTR-S	20 µL clear filtered tips, sterile	96 tips/rack, 24 racks/case
	TTF-50-C-HTR-S	25 µL clear filtered tips, sterile	96 tips/rack, 24 racks/case
	TT-50-C-HTR	50 µL clear non-filtered tips	96 tips/rack, 24 racks/case
	TT-50-C-HTR-36	50 µL clear non-filtered tips	96 tips/rack, 2 racks/pack, 36 packs/case
	TTF-200-C-HTR-S	175 µL clear filtered tips, sterile	96 tips/rack, 24 racks/case
	TT-200-C-HTR	200 µL clear non-filtered tips	96 tips/rack, 24 racks/case
	TT-200-C-HTR-S	200 µL clear non-filtered tips, sterile	96 tips/rack, 24 racks/case
	TT-1000-C-HTR	1000 µL clear non-filtered tips	96 tips/rack, 16 racks/case
	TT-1000-C-HTR-S	1000 µL clear non-filtered tips, sterile	96 tips/rack, 16 racks/case
	TTF-1000-C-HTR-S	1000 µL clear filtered tips, sterile	96 tips/rack, 16 racks/case

### Axygen® 96-well Format Tips, Black Conductive, Liquid Sensing

Compatible with Tecan LiHa Head

- ▶ ANSI/SLAS hanging tip, polypropylene rack



Cat. No.	Description	Packaging
TT-20-CBK-HTR	20 µL black non-filtered tips	96 tips/rack, 24 racks/case
TTF-20-CBK-HTR-S	20 µL black filtered tips, sterile	96 tips/rack, 24 racks/case
TTF-50-CBK-HTR-S	25 µL black filtered tips, sterile	96 tips/rack, 24 racks/case
TTF-50-CBK-HTR-S36	25 µL black filtered tips, sterile (economy package)	96 tips/rack, 2 racks/pack, 36 packs/case
TT-50-CBK-HTR	50 µL black non-filtered tips	96 tips/rack, 24 racks/case
TT-50-CBK-HTR-36	50 µL black non-filtered tips (economy package)	96 tips/rack, 2 racks/pack, 36 packs/case
TT-50-CBK-HTR-S	50 µL black non-filtered tips, sterile	96 tips/rack, 24 racks/case
TTF-200-CBK-HTR-S	175 µL black filtered tips, sterile	96 tips/rack, 24 racks/case
TT-200-CBK-HTR	200 µL black non-filtered tips	96 tips/rack, 24 racks/case
TT-200-CBK-HTR-S	200 µL black non-filtered tips, sterile	96 tips/rack, 24 racks/case
TT-1000-CBK-HTR	1000 µL black non-filtered tips	96 tips/rack, 16 racks/case
TT-1000-CBK-HTR-24	1000 µL black non-filtered tips (economy package)	96 tips/rack, 24 racks/case
TT-1000-CBK-HTR-S	1000 µL black non-filtered tips, sterile	96 tips/rack, 16 racks/case
TTF-1000-CBK-HTR-S	1000 µL black filtered tips, sterile	96 tips/rack, 16 racks/case



### Axygen Single-well Automation Reservoirs (ANSI/SLAS Footprint)

Cat. No.	Description	Qty/Cs
RES-SW1-LP	Reagent reservoir with 1-channel trough, low profile	25
RES-SW1-LP-SI	Reagent reservoir with 1-channel trough, low profile, sterile, individually wrapped	25
RES-SW8-HP	Reagent reservoir with 8-channel trough, high profile	25
RES-SW8-HP-S	Reagent reservoir with 8-channel trough, high profile, sterile	25
RES-SW8-HP-SI	Reagent reservoir with 8-channel trough, high profile, sterile, individually wrapped	25
RES-SW12-HP	Reagent reservoir with 12-channel trough, high profile	25
RES-SW12-HP-SI	Reagent reservoir with 12-channel trough, high profile, sterile, individually wrapped	25
RES-SW96-LP	Reagent reservoir with 96-channel trough, low profile	25
RES-SW96-LP-SI*	Reagent reservoir with 96-channel trough, low profile, sterile, individually wrapped	25
RES-SW96-HP-M*	Reagent reservoir with 96-channel trough, high profile	25
RES-SW96-HP-SI-M*	Reagent reservoir with 96-channel trough, high profile, sterile, individually wrapped	25
RES-SW384-LP	Reagent reservoir with 384-channel trough, low profile	25
RES-SW384-LP-SI	Reagent reservoir with 384-channel trough, low profile, sterile, individually wrapped	25
RES-SW384-FX	Reagent reservoir with 384-channel trough, medium profile	25
RES-SW384-FX-SI	Reagent reservoir with 384-channel trough, medium, sterile, individually wrapped	25
RES-SW384-HP-SI	Reagent reservoir with 384-channel trough, high profile, sterile, individually wrapped	25

\*Made to order. Please contact your local Corning Account Manager for more details.

### Axygen Multiple Well Automation Reservoirs (ANSI/SLAS Footprint)

Cat. No.	Description	Qty/Cs
RES-MW4-HP	Reagent reservoir with 4-channel trough, high profile	25
RES-MW8-LP	Reagent reservoir with 8-channel trough, low profile	25
RES-MW8-LP-SI	Reagent reservoir with 8-channel trough, low profile, sterile, individually wrapped	25
RES-MW8-HP	Reagent reservoir with 8-channel trough, high profile	25
RES-MW8-HP-SI	Reagent reservoir with 8-channel trough, high profile, sterile, individually wrapped	25
RES-MW12-LP	Reagent reservoir with 12-channel trough, low profile	25
RES-MW12-LP-SI	Reagent reservoir with 12-channel trough, low profile, sterile, individually wrapped	25
RES-MW12-HP	Reagent reservoir with 12-channel trough, high profile	25
RES-MW12-HP-SI	Reagent reservoir with 12-channel trough, high profile, sterile, individually wrapped	25