

# CaliMat

## New pH Buffer Solutions

pH value (20 °C)	Volume	Order no.
2.00 ± 0.02	250 ml	CS-P0200/250
4.00 ± 0.02	250 ml 1000 ml 3000 ml	CS-P0400/250 CS-P0400/1000 CS-P0400/3000
7.00 ± 0.02	250 ml 1000 ml 3000 ml	CS-P0700/250 CS-P0700/1000 CS-P0700/3000
9.00 ± 0.02	250 ml 1000 ml 3000 ml	CS-P0900/250 CS-P0900/1000 CS-P0900/3000
12.00 ± 0.05	250 ml	CS-P1200/250
<b>Buffer sets</b>		
Set 4.00	3 x 250 ml	CS-PSET4
Set 7.00	3 x 250 ml	CS-PSET7
Set 9.00	3 x 250 ml	CS-PSET9
Set 4.00, 7.00, 9.00	je 250 ml	CS-PSET479



- Ready to use pH buffer solutions with a wide range of values pH 2.00 / 4.00 / 7.00 / 9.00 / 12.00
- Traceable to SRM / NIST
- Reasonable
- Sure through colorcoding
- For use in industrial process applications as well as in laboratories
- Convenient and safe handling thanks to integrated calibration beaker
- Ideal filling quantity of 3000 ml for automatic calibration systems with Unical 9000
- Supplied with quality control certificate